

Performance Monitor

Quick Start Guide



Works with
BodiMetrics Mobile App
FOR APPLE AND ANDROID

Questions? Contact us.


Toll Free: (844) 744-8800

Email: info@bodimetrics.com



BodiMetrics LLC
1601 N.Sepulveda Blvd., Suite 839
Manhattan Beach, California 90266
www.BodiMetrics.com

Welcome to BodiMetrics Performance Monitor!



WARNING REMINDERS



This product is not suitable for people wearing a pacemaker or any other internal medical device.



This product is not a medical device and should not be used to diagnose or treat any medical condition.



Do not use device in direct sunlight or in very bright fluorescent lighted areas; no or inaccurate results may occur.

General Navigation

Your device has two main menu screens. Use arrows on right and left of the screen to move between menu screens.



Settings



- Press the **Home** button.
- Navigate to second screen by pressing right arrow.
- Press the the **Settings** icon.
- Press the **Date and Time** menu to enter date and time of day. Press arrow at bottom of screen to complete your entry and return to the Settings menu.
- Press the **User Management** menu, then select to add a new profile. Enter information then press to complete.
- Press the **Voice Guide** option to turn ON/OFF voice prompts. Press the **Volume** chart 1 to 3 times to raise or lower volume of voice guide.
- On your next visit to your physician's office, or if you have a blood pressure cuff available, you may calibrate your systolic blood pressure as follows.

Press the **Settings** icon. Press down arrow and then select **SBP Calibration** from the menu. Follow on screen and voice instructions. You need to perform 2 calibrations and the device will inform you of successful completion. Once you calibrate your blood pressure, you will be able to track your systolic blood pressure. **Note:** please recalibrate your blood pressure every 90 days.



Heart Activity




Press the **Heart Activity** icon on your device screen.

- Before beginning, sit down, relax, hold device gently as shown and rest your arms on a table.
- Stay relaxed during measurement as exerting too much pressure or actively moving may affect conductivity.
- Dry skin or thick calluses may also affect conductivity. Moisten fingers and/or palm to improve conductivity when necessary.



Pulse Oximeter




Press the **Pulse Oximeter** icon on your device screen.

- Gently insert your right index finger under the plastic flap onto the internal SpO₂ sensor and relax.
- Do not press down too hard on sensor as it may impact your normal blood flow and cause inaccurate results.



Body Check

Press the **Body Check** icon on your device screen.



- Before beginning, sit down, relax, hold device gently as shown and rest your arms on a table.
- Stay relaxed during measurement as exerting too much pressure or actively moving may affect conductivity.
- Dry skin or thick calluses may also affect conductivity. Moisten fingers and/or palm to improve conductivity when necessary.



Temperature

Press the **Temperature** icon on your device screen.





- Press the **Body** icon on your device screen.
- Put the **temperature sensor** on temple as shown.
- Press and release the **home** button, then scan toward the forehead.
- The measurement is finished when you hear the **loud** beep.



Relax Me

Press the **Relax Me** icon on your device screen.



- Press **User** name or add a new one.
- Before beginning, sit down, relax, hold device gently as shown and rest your arms on a table.
- Stay relaxed during measurement as exerting too much pressure or actively moving may affect conductivity.
- Follow the breathing pattern as shown above, inhaling as the indicator moves upward and exhaling as it moves downward.
- Press the **Home** button to stop the measurement at any time. You will receive your Relaxation Index upon stopping the measurement.



Pedometer

Press the **Pedometer** icon on your device screen.

- Press the **Flag** icon to set a "goal" distance for your steps and fat burn
- Press the **Home** button to stop and calculate.
- Press the **Reports** icon on your device screen.



Reports

Press the **Reports** icon on your device screen.

- Press the individual report icon that you want to review.
- Select the **User** you want to track and the dated report.
- Follow on screen instructions.



Reminders

Press the **Reminders** icon on your device screen.

- Press **User** name or add a new one under **Settings** for new users.
- Press then enter a task or event. Select day and month of year. **Note:** You can add a new task at any time.



Bluetooth

Press the **Bluetooth** icon on your device screen.

- To communicate between your device and a smartphone you will need to open your BodiMetrics mobile app.
- **Note:** Turn off your Bluetooth when not in use, as it will extend your device's battery life.